
# **X-MIND<sup>®</sup> AC**

Intra-oral x-ray system

## **INSTALLATION & MAINTENANCE MANUAL**

**CE**  
**0434**

**This manual should always be kept near the installation**

*03/2001 Edition*

## SUMMARY

SUMMARY	page	2
INTRODUCTION	page	3
PRELIMINARY INFORMATION	page	4
INFORMATION FOR THE FITTER	page	5
1. "X-Mind® AC" RADIOGRAPHIC SYSTEM	page	6
2. IDENTIFICATION TAGS	page	7
3. OVERALL DIMENSIONS	page	8
4. INSTALLATION SPECIFICATIONS	page	12
5. INSTALLATION	page	15
5.1 UNPACKING	page	15
5.2 ASSEMBLING THE WALL PLATE	page	16
5.3 ASSEMBLING THE BRACKET	page	18
5.4 ASSEMBLING THE PANTOGRAPH TYPE ARM	page	20
5.5 CONNECTING TO THE FEEDING TERMINAL BOARD	page	21
5.6 ASSEMBLING THE RADIOGRAPHIC UNIT	page	24
5.7 BALANCING THE PANTOGRAPH TYPE ARM	page	27
5.8 ASSEMBLING THE SWITCHBOARD	page	30
5.9 SWITCHBOARD ELECTRICAL CONNECTION	page	32
6. CONTROL PANEL	page	34
7. SYSTEM CONFIGURATION	page	36
8. CHANGING THE CONFIGURATION	page	37
9. START UP	page	38
10. CHECKING THE INSTALLATION	page	39
11. CHECKING THE EXPOSURE FACTORS	page	41
12. DIAGNOSTIC	page	44
13. FAULT MESSAGES	page	45
14. REPLACING FUSES	page	46
15. CLEANING THE OUTER SURFACES	page	47
16. MAINTENANCE	page	48
17. REPAIRS	page	48
18. ATTACHMENTS	page	50

## INTRODUCTION

The radiographic system described in this manual is a "**wall installation**".

SATELEC® S.A.S. reserves the right to modify its products and manual without notice.

SATELEC® S.A.S. shall not be liable for any incorrect use of the information contained in this manual.

Any copies, even partial, of this manual are permitted solely for in-house use.

## PRELIMINARY INFORMATION

Before using the "X-Mind® AC" radiographic system, we recommend you carefully read and follow the instructions contained herein, to get the maximum out of the equipment.

Always pay close attention to the **CAUTION**, **WARNING**, and **PLEASE NOTE** messages when operating the system.

### LEGEND

#### **CAUTION**

The word **CAUTION** identifies possible incidents, which could endanger **the operator's personal safety or cause personal injuries**.

#### **WARNING**

The word **WARNING** identifies those incidents, likely to affect **the radiographic system's performance**.

#### **PLEASE NOTE**

The words <**PLEASE NOTE**> are used to highlight particular points to facilitate maintenance or make important information clearer.

## INFORMATION FOR THE FITTER

### **CAUTION**

***The fitter is responsible for the installation, the system safety and operation.***

To ensure that the "X-Mind® AC" radiographic system is correctly and safely installed, it is advisable to:

- ✓ Check that the voltage shown on the rating plate matches the line voltage
- ✓ Install the "X-Mind® AC" radiographic system following to the procedures described in this manual.
- ✓ Provide the user with all information regarding the use of the "X-Mind® AC" radiographic system according to what is stated in the manual
- ✓ Provide a "Declaration of Conformity" to certify the work done
- ✓ Return to "SATELEC S.A.S." the warranty certificate duly filled in:  
Failure to do so will render the warranty invalid


## 1. "X-Mind® AC" RADIOGRAPHIC SYSTEM

The "X-Mind® AC" radiographic system (Fig. 1) consists of:

- ① Tube head
- ② Spacer cone
- ③ Pantograph type arm
- ④ Switchboard
- ⑤ Wall plate
- ⑥ Bracket


### OPTIONAL

- short 8" (20cm) cone
- cone with a rectangular section sized 44x35mm
- second "CONTROL BUTTON" with extendable cable
- RX signalling lamp for external use


## 2. IDENTIFICATION TAGS

The identification tags on the tube head, on the switchboard and on the cone indicate the model number, the serial number, the manufacturing date and the symbols of the main technical characteristics.

ID TAG OF THE  
TUBE HEAD

ID TAG OF THE  
SWITCHBOARD

ID TAG OF THE LONG 12" CONE

ID TAG OF THE SHORT 8" CONE

ID TAG OF THE RECTANGULAR CONE

GRADUATED SCALE TAG ON TUBE HEAD

### Pictograms used

This symbol guarantees that the radiographic system complies with the regulations of the European Directive EEC 93/42 on Medical Devices

The degree of protection against direct and indirect electrical contacts is B type

Refer to the instructions in the manual

**N**

"NEUTRAL" power supply wire

**L**

"PHASE" power supply wire


Earth cable

### 3. OVERALL DIMENSIONS


Fig. 1A, 1B, 1C give the overall dimensions of the possible supply conditions:

☛ Short bracket (optional)	length 41cm	- 16,2"
☛ Standard bracket	length 82,5cm	- 32,5"
☛ Long bracket (optional)	length 110cm	- 43,5"


SHORT BRACKET 41cm

*Fig. 1A*

STANDARD BRACKET 82,5cm

*Fig. 1B*


LONG BRACKET 110cm

*Fig. 1C*

Fig. 2, 3A, 3B, 3C and 4 show the typical dimensions of the radiographic system:

STANDARD BRACKET 82,5cm

Fig. 2

STANDARD BRACKET 82,5cm

Fig. 3A



Fig. 3B  
SHORT BRACKET 41cm


Fig. 3C  
LONG BRACKET 110cm


Fig. 4

## 4. INSTALLATION SPECIFICATIONS

### **WARNING**

*Prior to installing the radiographic system the Office Owner must ensure that: **the environment, the electrical system and the power supply meet the requirements needed.** If this is not the case, he shall make the necessary adjustments.*

### **1. ENVIRONMENT REQUIREMENTS**

- a. The installation environment must be of suitable width.  
With the sizes and overall dimensions provided (see §3 "OVERALL DIMENSIONS"), check that there are no obstacles when positioning the radiographic system
- b. The environment must not be exposed to risks of explosion and must not be pressurized
- c. During operation, the ambient temperature range must be between +5°C and + 40°C
- d. The storage temperature range must be between - 15°C and + 50°C
- e. The relative humidity must range between 25% and 75%

### **2. SUPPORTING WALL REQUIREMENTS**

- a. The radiographic system supporting wall must be able to support a load of 200Kg at each fixing point
- b. The nature and consistency of the wall must be checked. If required, consult a construction expert
- c. A buried counter plate or sandwich type system shall be fitted to walls of uncertain consistency

### **3. ELECTRICAL SYSTEM REQUIREMENTS**

- a. The electric system must comply with the regulations in force
- b. The electric system must be able to supply the power and voltage required in the Manufacturer's rating plate of the radiographic system as per Chart 1

**CHART 1**

<b>Radiographic system input voltage</b>	<b>220 V -5% +10%</b>	<b>230 V -5% +10%</b>	<b>240 V -5% +10%</b>	<b>115 V -5% +10%</b>
Nominal line voltage	220V	230V	240V	115V
Minimum line voltage	209 Vac	218.5 Vac	228 Vac	109 Vac
Maximum line voltage	242 Vac	253 Vac	264 Vac	126,5 Vac
Frequency	50/60 Hz	50/60 Hz	50/60 Hz	50/60 Hz
Power absorbed during exposure (max. time = 3,9sec.)	800 VA	800 VA	800 VA	800 VA

#### 4. ELECTRICAL LINE REQUIREMENTS

- a. The electrical line must be "single-phase" type
- b. It is essential a 10A – 250V, magneto thermal differential switch be fitted at the start of the radiographic system, with differential protection  $I_n \leq 30\text{mA}$   
(see §18 "ATTACHMENTS")
- c. The power conductors of the switchboard and the connection conductors to the radiographic unit must be two-pole + earth with suitable section to the length of the power supply line, as per Chart 2

**CHART 2**

Power supply voltage	Minimum conductor section	Maximum line length
209 to 264V	1.5 mm <sup>2</sup>	10 metres
	2.5 mm <sup>2</sup>	20 metres
109 to 126.5V	1.5 mm <sup>2</sup>	5 metres
	2.5 mm <sup>2</sup>	10 metres

#### PLEASE NOTE

For longer lines, the conductor section must be increased proportionately.

- d. The cables connecting the switchboard and the signal lamps located outside the office must be two-pole type, of section  $\geq 0.5 \text{ mm}^2$
- e. The electrical line characteristics must comply with Chart 3

**CHART 3**

Line voltage	220 V -5% +10%	230 V -5% +10%	240 V -5% +10%	115 V -5% +10%
Maximum line voltage drop	3%	3%	3%	3%
Maximum line apparent resistance	0.5 ohm	0.5 ohm	0.5 ohm	0.2 ohm

## 5. ELECTRICAL CONNECTIONS

### **WARNING**

*Prior to installing the radiographic system, we recommend that all the electrical connections be prepared.*

### **Switchboard electrical connections**

### **CAUTION**

***In accordance with the relevant standard, the switchboard must be installed in a position allowing the operator to permanently control the radiographic exposure.***

The switchboard installation wall should be fitted with suitable runs for the following electrical cables, as per the installation electrical diagram (see §18 "ATTACHMENTS"):

- a. Switchboard electrical cables  
(see §5.9 "SWITCHBOARD ELECTRICAL CONNECTIONS")
- b. Cables connecting the switchboard and the radiographic unit  
(see §5.9 "SWITCHBOARD ELECTRICAL CONNECTIONS")
- c. Cables connecting the switchboard and the signal lamps outside the office  
(IF PROVIDED)

### **Radiographic unit electrical connections**

The wall plate installation wall should be fitted with a suitable run for the cable connecting the switchboard and the radiographic unit  
(see §5.5 "CONNECTION TO THE POWER SUPPLY TERMINAL BOARD")

## 5. INSTALLATION

### CAUTION

***Only professionally trained technicians able to certify their work with a "Declaration of Conformity" can fit the "X-Mind® AC" radiographic system***


### CAUTION

***Prior to commencing installation, ensure that all the requirements have been met.***

### 5.1 UNPACKING THE EQUIPMENT

The components of the "X-Mind® AC" radiographic system come packed inside a cardboard box, as shown in the following sketch (Fig. 5):

- D Documents: Instructions Manual, Installation and Maintenance Manual, Warranty
- ① Radiographic unit packaging
- ② Switchboard packaging+ External Lamp (IF PROVIDED) packaging
- ③ Pantograph type arm packaging
- ④ Bracket type wall plate packaging

### PLEASE NOTE

Check all components prior to installation

### PLEASE NOTE

The cardboard box and the polystyrene foam can be completely recycled and can be sent for disposal to authorized recycling companies.

### PLEASE NOTE

We recommend keeping the original packaging in case goods need to be returned for repair.

## 5.2 FITTING THE WALL PLATE

### **WARNING**

*DO NOT use plastic or rubber anchor screws to fit the wall plate.  
Use metal anchor screws Ø12 (NOT included in the supply) when fitting to cement walls, or walls built with solid or hollow bricks.*

### **ASSEMBLY INSTRUCTIONS (Fig. 6)**

1. Remove the wall plate from the packaging (see Fig. 5) and take the drilling template ①
2. Position the drilling template ① on the radiographic system installation wall, at the required height (130cm from the base is the recommended height)
3. Fix the template using adhesive tape
4. Using a plumb line, check that the holes are perpendicular and aligned with the floor.
5. Mark out the plate fixing holes
6. Mark out the holes for the electrical cables connecting the switchboard to the radiographic unit

### **PLEASE NOTE**

To prevent any flaking in the white coat and to control the centre distances between the holes, it is advisable to start drilling with a tip Ø7, and increase this measure gradually.

7. Drill the plate fixing holes
8. If required, drill the holes for the electrical cables connecting the switchboard to the radiographic unit
9. Remove the template ① and insert anchor screws suiting the wall characteristics ②
10. Unscrew the screw ④ and remove the plug ③ from the wall plate
11. Withdraw the sliding cover ⑤
12. Set the plate ⑥ to the wall and insert the screws ⑦ with the relevant washers, then tighten alternately
13. Check that the plate is fixed firmly to the wall

### **PLEASE NOTE**

If the wall is not completely level, put a suitable shim between the wall and the plate, to prevent any deformation.


Fig. 6


### 5.3 ASSEMBLING THE BRACKET

#### PLEASE NOTE

The brackets come in the following lengths: 41cm - 82.5cm - 110cm  
The 82.5cm and 110cm brackets are fitted with a stop key ① (Fig. 7A and 7B) to prevent the electrical cable from twisting.


#### PLEASE NOTE

Generally, the stop key is fitted so that the rest position of the equipment is to the right of a possible observer standing in front of the wall plate (Fig. 7A).  
If the equipment comes to rest to the left, the stop key must be rotated 180° (Fig. 7B).

After fitting the plate to the wall, assemble the bracket using the following procedure (Fig. 8):

1. Remove the bracket from the packaging (*see Fig. 5*)
2. Insert the bracket pin ② into the wall plate ① (upwards)
3. Insert the supporting rest ③

*Fig. 8*

**PLEASE NOTE**

Prevent any foreign bodies (earth, dust, cement, etc.) from settling on the pin seat. The pin should slide freely into its seat. If required, clean thoroughly and lubricate with "Molikote D" grease.

**PLEASE NOTE**

Carefully check the bracket is parallel on the floor using a spirit level.

## 5.4 ASSEMBLING THE PANTOGRAPH TYPE ARM

### ASSEMBLY INSTRUCTIONS (Fig. 9)

1. Remove the pantograph type arm from the packaging (*see Fig. 5*)
2. Remove the bracket plug ① by unscrewing the fixing screw ②
3. Take off the bracket guard slat ③
4. Insert the pantograph group cable ⑤ into the washer ④ and then the pantograph pin ⑥


Fig. 9

5. If required, clean the pin and the washer and lubricate with "Molikote D" grease
6. Insert the electric cable into the bracket housing ⑦
7. Reposition the guard slab
8. Insert the cable into the bracket. Push it until it reaches the pin outlet near the supply terminal board

## 5.5 CONNECTION TO THE FEEDING TERMINAL BOARD

### **CAUTION**

*For electrical safety, it is essential that the earth wires be correctly connected.*

### **WARNING**

*The electrical cable length is suitable for 82.5cm long brackets.  
DOT NOT tamper with the terminal arrangement when 41cm long brackets are used.  
The excess cable must be housed in the housing provided.  
The 110cm long brackets come with a special extension cable.*


### **WARNING**

*In the case of PHASE – PHASE type electrical power supply, the higher value must be taken using a tester and entered in L.*

## OPERATING INSTRUCTIONS (Fig. 10)

1. Remove the terminal board cover ① by unscrewing the fixing screw ②

*Fig. 10*

2. Complete the electric connection as shown in Fig. 11


Fig. 11

#### TERMINAL BOARD CONNECTION DIAGRAM

①	BROWN	L (line)
②	YELLOW GREEN	T (earth)
③	BLUE	N (neutral)
④	NOT CONNECTED	
⑤	NOT CONNECTED	
⑥	YELLOW GREEN	T (earth)
⑦	YELLOW GREEN	T (earth)

3. Connect the pantograph cable shield to the earth potential ⑥
4. Connect the wall plate to the earth potential ⑦
5. Clamp the cables with the cable clamps provided
6. Refit the terminal board insulating cover

## 5.6 ASSEMBLING THE RADIOGRAPHIC UNIT

### ASSEMBLY INSTRUCTIONS

1. Remove the radiographic unit from the packaging (*see Fig. 5*)
2. Check that all the rating data match the power supply voltage
3. Remove both guards from the pantograph type arm by loosening the relevant screws (Fig. 12)

*Fig. 12*

4. Using a tip, work on the front coupling device (Fig. 13A)
5. Remove both guards (Fig. 13B)

*Fig. 13A*


*Fig. 13B*

6. Insert the radiographic unit pin into the pantograph head (Fig. 14A) and insert the supporting rest (Fig. 14B)

*Fig. 14A*

*Fig. 14B*

7. Check that during insertion, the pin of the anti-spinning device fits correctly into the seat located on the pantograph head (Fig. 15)

*Fig. 15*

8. Couple the pantograph and mono-bloc connectors and fit them into their seats (Fig. 16)

*Fig. 16*

## 5.7 BALANCING THE PANTOGRAPH TYPE ARM

### CAUTION

*The pantograph type arm must only be adjusted when the radiographic unit is assembled.*

### WARNING

*To prevent damage to the internal mechanism, the adjustment key must not be left in place when adjustments and balancing tests are being carried out.*

### WARNING

*Do not lose the key provided.*

### PLEASE NOTE

To reach the adjustment screw X, set the arm "A" to the vertical position.

To reach the adjustment screw Y, set the arm "B" to the horizontal position.

The adjustment key provided can only be inserted under the above conditions (Fig. 17)


Fig. 17

Procedure to be followed to correctly balance the pantograph type arm (*see Fig. 17*):

## **1. BALANCING THE ARM "A"**

### **PLEASE NOTE**

The pantograph type arm is supplied with arm "A" already preloaded.  
The Arm "B" comes unloaded for safety reasons.

## **2. BALANCING THE ARM "B"**

For correct balancing, proceed as follows:

1. Arm "A" vertical
2. Arm "B" horizontal
3. Insert the adjustment key into "Y"
4. Load the spring turning it 22 times
5. **Remove the key**

## **3. CHECKING THE BALANCING**

To check the balancing proceed as follows:

set the arm "B" to the different positions

if it does not keep the correct position:

1. Set the arm "B" back to the horizontal position
2. Insert the adjustment key into "Y"
3. Rotate the adjustment screw by a half turn:
  - clockwise if it tends to come down
  - anticlockwise if it tends to go up
4. **Remove the key**

### **PLEASE NOTE**

Repeat the test and adjustment until the arm "A" is steady and stable in all positions, even when the arm "B" is completely extended.

## **4. READJUSTMENT OF ARM "A"**

If it does not stay in the required position:


1. Set the arm "A" to the vertical position again
2. Insert the adjustment key into "X"
3. Turn the adjustment key by a half turn
4. clockwise if it tends to come down
5. anticlockwise if it tends to come down
6. **Remove the key**

### **PLEASE NOTE**

Repeat the test and adjustment until the arm "A" is steady and stable in all positions, even when the arm "B" is completely extended.


Once balancing is complete:

1. Insert the movable finish (A) between the pantograph guard and metal frame (Fig. 18)

*Fig. 18*

2. Insert the pins ① of one guard ② into the relevant seats; then position it and check that the finish ③ is coupled to the guards (Fig. 19)

*Fig. 19*

## 5.8 ASSEMBLING "X-Mind® AC" SWITCHBOARD

### CAUTION

**Check that the cable runs are arranged in the switchboard installation wall and check that the power supply meets the installation specifications (see §.4 "INSTALLATION SPECIFICATIONS").**


### CAUTION

**Check that the rating data match the power supply voltage.**

### ASSEMBLY INSTRUCTIONS (Fig. 20)

1. Remove the switchboard from the packaging (see Fig. 5)
2. Take out the drilling template ①
3. Mark the fixing holes on the wall using the drilling template
4. Drill using a Ø3mm tip, then drill again with a Ø6mm tip to prevent any flaking on the white coat
5. Insert the three screw anchors provided into the holes ②
6. Open the switchboard by unscrewing the three screws ③

Fig. 20

7. Remove the 26-pole connector from its seat to release both switchboard guards (Fig. 21)


Fig. 21

8. Insert the electric feeding cables into the hole (Fig. 22)


Fig. 22

9. Insert the connection cables coming from the radiographic units into the rectangular slot
10. Insert the cables of the external signal lamps (if installed) (Fig. 22)
11. Set the switchboard base to the wall, matching the three anchor screws with the holes
12. Screw the screws (*see Fig. 20 ④*) with the relevant washers (*see Fig. 20 ⑤*) into the screw anchors

**PLEASE NOTE**

If the wall is uneven, set in some shims to prevent the switchboard from buckling.

**PLEASE NOTE**

Do not soil the switchboard with dust or rubble resulting from drilling.

## 5.9 SWITCHBOARD ELECTRICAL CONNECTION

### CAUTION

***Before carrying out any connections, ensure that the power supply is isolated.***

### CAUTION

***If the equipment is mounted on metal walls, the latter must be connected to the ground circuit.***

### WARNING

***Always comply with the polarity PHASE – NEUTRAL when making connections.***

### WARNING

***When stripping the cables, pay attention to the small copper wires that may fall on the printed circuit and cause short circuits or malfunctioning.***

### OPERATING INSTRUCTIONS (Fig. 23) (see §18 "ATTACHMENTS")


1. Connect the power supply cable to the terminal board ⑥
2. Insert the three mains cables into the rack
3. Fix them using the cable clamp ⑦
4. Connect the cables coming from the Radiographic unit 1 to the terminals XRAY1
5. Connect the YELLOW-GREEN earth cable to the "equipotential metal plate" ⑧
6. Connect the cables coming from the Radiographic unit 2 to the terminals XRAY2
7. Connect the YELLOW-GREEN earth cable to the "equipotential metal plate" ⑧
8. Clamp the cables in the cable clamp ⑨
9. Connect the RX signal lamps (if installed) (see the relevant "PLEASE NOTE")
10. Check the configuration on the dip-switches (see §7 "SYSTEM CONFIGURATION")
11. Reconnect the 26-pole connector
12. Close the switchboard with the three screws
13. Refit the plate guard (see Fig. 6 ⑤)
14. Reconnect the power supply

### PLEASE NOTE

To connect the RX warning lamps outside the office, see the instructions supplied with the relevant packaging.

The electronic card terminal boards have the following functions:

- 1 RADIOGRAPHIC UNIT No. 1 CONTROL BUTTON
- 2 RADIOGRAPHIC UNIT No. 2 CONTROL BUTTON (OPTIONAL)
- 3A RADIOGRAPHIC UNIT No. 1 SIGNAL LAMP OUTSIDE THE OFFICE
- 4A RADIOGRAPHIC UNIT No. 2 SIGNAL LAMP OUTSIDE THE OFFICE
- 3B NOT AVAILABLE
- 4B NOT AVAILABLE
- 5 RADIOGRAPHIC UNIT No. 1 POWER SUPPLY
- 6 RADIOGRAPHIC UNIT No. 2 POWER SUPPLY
- 7 SWITCHBOARD POWER SUPPLY

## 6. CONTROL PANEL


The control panel is the user's interface, where the operator can:

using the keys

- ➡ select the radiographic unit operating parameters
- ➡ change the exposure times displayed
- ➡ store the new customised exposure times


with the light indicators "green Led"

- ➡ check the switchboard operating conditions
- ➡ check the active function key
- ➡ check the selected radiographic unit parameter

## 7. SYSTEM CONFIGURATION


The "X-Mind® AC" radiographic system is factory configured for an operative "standard mode" which determines:

Installation of 2 "X-Mind® AC" tube heads	⇒	On the control panel by pressing the RX button  Led 1 lights up. By pressing the button again, Led 2 lights up.
A radiographic distance DSS = 31cm with long 12" cone	⇒	on the control panel Led 12"  is lit
with film type "D"	⇒	on the control panel led "D"  is lit
CONTROL BUTTON to perform the exposure	⇒	The switchboard houses a key with extendable cable

The configuration may be changed if:

<ul style="list-style-type: none"> <li>– type "E" films are used</li> <li>– a digital system is used</li> </ul>	⇒	by pressing the keys of the control panel <i>see USER'S MANUAL §6 "USE INSTRUCTIONS"</i>
<ul style="list-style-type: none"> <li>– the short 8" (20cm) cone is used</li> <li>– one single radiographic unit is used</li> <li>– 2 CONTROL BUTTON are used</li> </ul>	⇒	by changing the dip-switch position <u>ONLY THE FITTER MUST CARRY OUT THIS OPERATION</u>

The installed radiographic units have the following features:

They work with alternating current		the Led "AC" on the control panel must be lit
They work with a radiographic voltage equal to 70kVp		the Led "70 kVp" on the control panel must be lit
They work with a radiographic current equal to 8ma		the Led "8 mA" on the control panel must be lit

The above configuration depends on the position of n° 8 mini-switches (dip-switch) on the switchboard electronic card.

	ON	OFF	
1	■		radiographic unit No. 1
2		■	not available
3	■		radiographic unit No. 2
4		■	not available
5		■	control button No. 2
6	■		cone
7		■	not available
8		■	not available

## 8. CHANGING THE CONFIGURATION

### CAUTION

*The fitter must carry out this operation.*

### WARNING

*To make the changes operative, turn the switchboard off and then turn it on again.*

To change the configuration, the dip-switch positions in the switchboard must be changed:

DIP SWITCH	PARAMETER	ON	OFF
1	RADIOGRAPHIC UNIT No. 1	INSTALLED	NOT INSTALLED
2	TYPE OF RADIOGRAPHIC UNIT No. 1	NOT AVAILABLE	
3	RADIOGRAPHIC UNIT No. 2	INSTALLED	NOT INSTALLED
4	TYPE OF RADIOGRAPHIC UNIT No. 2	NOT AVAILABLE	
5	CONTROL BUTTON No. 2	INSTALLED	NOT INSTALLED
6	CONE	LONG (12") INSTALLED	SHORT (8") INSTALLED
7	NOT AVAILABLE	NOT AVAILABLE	
8	NOT AVAILABLE	NOT AVAILABLE	

### 1. CHANGING THE AMOUNT OF RADIOGRAPHIC UNITS INSTALLED

To change the amount of radiographic units installed, move the dip-switch No. 1 or No. 3

- If the radiographic unit is connected to the RX1 terminal board, set dip-switch No. 1 to the ON position, otherwise to the OFF position
- If the radiographic unit is connected to the RX2 terminal board, set dip-switch No. 3 to the ON position, otherwise to the OFF position

### 2. CHANGING THE CONE

To replace the 12" cone with an 8" cone, move the dip-switch No. 6 to the OFF position. Check that the led "8" lights up on the control panel.

#### PLEASE NOTE

After the changes, the set exposure times are automatically changed.

### 3. CHANGING THE AMOUNT OF THE REMOTE BUTTONS

To change the amount of control buttons, move the dip-switch No. 5 to the ON position.

#### PLEASE NOTE

After the changes, the relevant button controls each radiographic unit.

## 9. START UP

### CAUTION


**For further information, see the USER MANUAL.**

**TURN ON THE SWITCHBOARD** to power up to the radiographic system

- a. Set the "KEY SWITCH" to the "I" position (ON)

- b. Set the "MAIN SWITCH" located on the upper part of the switchboard to the "I" position (ON)

- ➡ the green light comes on indicating that the system is powered
- ➡ the Leds of the preset radiographic parameters automatically light up
- ➡ the exposure time is shown on the display

- c. **THE RADIOGRAPHIC SYSTEM IS NOW READY FOR USE**

### CAUTION

***If a fault is detected when the system is turned on, the anomaly is indicated as follows:***

- ***an intermittent beep sounds***
- ***the "MALFUNCTIONING INDICATOR" Led  comes on intermittently***  
***The error code (E ....) appears on the display***  
***(see §13 "ERROR MESSAGES")***
- ***All CONTROL PANEL functions are inhibited***

***In this happens, turn off the switchboard and then turn it back on.***

***If the fault persists, call "technical support".***

### PLEASE NOTE

The exposure time and radiographic parameters appearing on the display are the last that were set before the switchboard was turned off.

### PLEASE NOTE

If installed, outside the office, the RX signalling lamp, corresponding to the selected tube head comes on.

### PLEASE NOTE

If the switchboard remains inactive for a few minutes, it switches to the stand-by mode. Press any key of the CONTROL PANEL to restore it to the operative mode.

## 10. CHECKING THE INSTALLATION

### CAUTION

***When all connections are completed, the fitter must check the electrical safety and functions of the system.***

## OPERATING INSTRUCTIONS

### 1. Checking the configuration


Check on the control panel that all LEDs corresponding to the required configuration are lit; otherwise, change them.

(see §7 "SYSTEM CONFIGURATION" and §8 "CHANGING THE CONFIGURATION")

### 2. Checking the switchboard operation

- a. Check that the control panel is working correctly by selecting different exposure times
- b. Check the time on the display
- c. Check that, when changing the selected radiographic unit, the corresponding external RX signal lamps (if installed) come on

### 3. Checking the exposure

- a. Set an exposure time of 1sec.
- b. Use the "CONTROL BUTTON" on the switchboard
- c. Using the extendable cable of the "CONTROL BUTTON", keep a safety distance of at least 2 metres from the radiographic unit, in order to be able to constantly check the radiographic exposure
- d. Press the "X-RAY KEY"  and keep it pressed until the acoustic signal (beep) stops and the yellow "X-RAY OUTPUT SIGNAL"  Led switches off
- e. At the end of the exposure, the green Led "PAUSE INDICATOR"  the pause period.
- f. Check the indication of the actual exposure time on the display  
(see §11 4. "EXPOSURE TIME (sec)")

### 4. Checking the operation of the radiographic unit

Carry out several exposures on the radiographic units installed and check that:

- a. There are no faults
- b. The Led of the selected radiographic unit is lit
- c. The "CONTROL BUTTON" Led is lit throughout the entire duration of the acoustic signal

## 5. Checking the power absorbed by the radiographic system

To check the power absorbed by the radiographic system a tester must be used, in the mode ammeter in "AC"

- a. Connect the instrument to the power supply line  
(see §5.9 "SWITCHBOARD ELECTRICAL CONNECTION")
- b. Set an exposure time of approx. 3sec. on the switchboard
- c. Carry out an exposure and read the current value on the instrument

### PLEASE NOTE

The radiographic system complies with the requirements when

the power absorbed is  $\leq 6.3A$  for line voltages of 220V 230V 240V  
the power absorbed is  $\leq 8A$  for line voltage of 115V

Otherwise, check the electric system or call "Technical Support".

## 6. Checking the electrical system

To check the electrical system, a tester must be used, in the mode Voltmeter in "AC"

- a. Connect the instrument to the terminals L and N on the switchboard
- b. Measure the line voltage (load voltage  $V_0$ )
- c. Connect the instrument to the terminals L and N of the wall plate terminal board (see §5.5 "CONNECTION TO THE FEEDING TERMINAL BOARD")
- d. Set an exposure time of approx. 3 sec. on the switchboard
- e. Carry out the exposure
- f. Measure the line voltage during exposure (load voltage  $V_1$ )

### PLEASE NOTE

The electrical system complies with the requirements when the result of the formula

$$(V_0 - V_1) / V_0 \text{ is } \leq 0.03 \text{ (3\%)}$$

Otherwise, the electrical system must be adjusted (see §4 "INSTALLATION SPECIFICATIONS").

## 11. CHECKING THE EXPOSURE FACTORS

### OPERATING INSTRUCTIONS

#### 1. Radiographic voltage (Kvp)

The radiographic high voltage is measured using a calibrated "non invasive" instrument with an initial delay of 50msec.

The kVp value is factory checked with an instrument PMX-I D.

##### **Set technical factors**

**Line voltage = V nominal – 5% +10%**

**Max. voltage drop = 3%**

**KVp = 70**

**mA = 8**

**Set exposure time = 1sec. – with a delay-time of 0.05sec**

The radiographic VOLTAGE is 70kVp, with +/-10% tolerance

#### 2. Radiographic current (mA)

The radiographic current is measured by connecting a milliammeter in parallel with the 1 kohm resistance, installed inside the radiographic unit.

The mA value is factory checked with a digital multimeter.

##### **Set technical factors**

**Line voltage = V nominal – 5% +10%**

**Max. voltage drop = 3%**

**KVp = 70**

**mA = 8**

**Set exposure time = 3.2sec.**

The radiographic CURRENT is = 8mA +/-10%

#### 3. Dose (mGy)

The dose in air is measured with a "non invasive" instrument, by positioning the detector at a spot-skin distance = 310mm

The mGy are factory checked with the instrument SOLIDOSE.

##### **Set technical factors**

**Line voltage = V nominal – 5% +10%**

**Max. voltage drop = 3%**

**KVp = 70**

**mA = 8**

**Set exposure time = 1sec.**

DOSE in air = 3.3 mGy/s +/-30%

#### 4. Exposure time (sec)

The exposure time is measured with a "non invasive" instrument.

##### PLEASE NOTE

The "X-Mind® AC" switchboard features "self-compensating" technology. The automatic correction or compensation of the exposure time gets round the effects of the line voltage fluctuations on the X-ray film quality. The radiographic over/under exposure is eliminated. The automatic correction works as follows:

- A voltmeter located inside the timer constantly reads the mains voltage, while the operator selects the required exposure time.
- The last measure collected before the exposure is sent to the internal microprocessor that corrects the exposure time set by the operator by applying a special algorithm.
- At the end of the exposure, during the pause time, the time corrected value is shown on the display of the switchboard.

THE PAUSE TIME IS APPROX. 32 TIMES THE SET TIME.

**The correct time is the actual duration of the exposure**

The accuracy of the "actual time" is factory checked with the instrument PMX-I.

##### Set technical factors

**Line voltage = V nominal**

**Max. voltage drop = 3%**

**KVp = 70**

**mA = 8**

**Set exposure time = 0.2sec**

The EXPOSURE TIME measured is = 0,2sec +/-10% or 20msec

#### 5. Checking the exposure time

##### **WARNING**

*Some differences between the set exposure time and the measured exposure time may be detected during periodic quality checks.*

*The differences are due to the automatic time compensation that comes into play when the power voltage differs from the nominal one.*

The switchboard comes with two types of accuracy measurement for the exposure time:

##### **A. BEFORE ACQUISITION**

(It is the exposure time set on the switchboard by the operator)

The max. error is +100% - 100% when the voltage fluctuation is -5% + 10%.

The max. error is  $\pm 10\%$  or  $\pm 20\text{msec}$  at nominal voltage.

##### **B. AFTER EXPOSURE**

(It is the exposure duration. It is calculated by the microprocessor according to the line voltage, and it is shown on the display)

The max. error is  $\pm 10\%$  or  $\pm 20\text{msec}$ . (1 pulse)

## OPERATING INSTRUCTIONS

To check the exposure times (sec) it is advisable to use calibrated "non invasive" instruments, i.e.:

NERO by Victoreen

PMX by RTI

TRIAD QA by Keithley

- a. Put the "X-ray detector" at the cone outlet, centred on the primary beam
- b. Turn on the switchboard
- c. Set one exposure time on the switchboard
- d. Carry out the exposure
- e. At the end of the exposure, when the functional pause period starts, read the exposure **actual time** on the display
- f. Use the chart below To check the accuracy of the exposure time,

### CONFORMITY CHECK CHART

At the end of the exposure

1° read the time value  
on the display

= Actual exposure time

2° read the measure  
on the instrument

= Emission time measure

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4° divide 3° by 2°

= Absolute error

5° multiply 4° by 100

= Percentage error (%)

The accuracy is checked when the  
% error is

by adopting the most favourable  
condition

### PLEASE NOTE

The exposure time complies with the requirements when:

$$((\text{Measured time} - \text{Actual time}) / \text{Measured time}) * 100 \leq \pm 10\% \text{ or } \pm 20 \text{ msec}$$

Otherwise, please contact the "Technical Support".

## 12. DIAGNOSTIC

The "X-Mind® AC" radiographic system allows the fitter to set or display some functional parameters of the switchboard:

1. To set the parameters, the fitter must:

- a. Turn off the switchboard
- b. Simultaneously press and hold down the keys  
(45) LOWER JAW PREMOLAR + (43) LOWER JAW CANINE
- c. Turn on the switchboard
- d. When the message "inst" is displayed, the installer can set the minimum exposure time:
  1. select the tube head
  2. turn off the switchboard
  3. turn on the switchboard, by holding down the  key
  4. the message "SEC" is displayed for approx 1sec
  5. the present value of the lower set limit is displayed
  6. to change the value, press the keys  and 
  7. to confirm, press the  key
- e. To exit this mode, turn the switchboard off and then on again

2. To visualise the parameters, the fitter must:

- a. Press simultaneously and hold down the keys  
(17) UPPER JAW MOLAR + (47) LOWER JAW MOLAR
- b. Press the key associated with the parameter to be viewed

KEY		DISPLAYED PARAMETER	Example	U.M.
	<b>BITE-WING ANT</b>	RADIOGRAPHIC SYSTEM VOLTAGE	220	Volt
	<b>BITE-WING POST</b>	LINE VOLTAGE	227	Volt
	<b>LOWER INCISORS</b>	MAXIMUM LINE VOLTAGE VALUE DETECTED	238	Volt
	<b>LOWER INCISORS</b>	MINIMUM LINE VOLTAGE VALUE DETECTED	215	Volt
	<b>OCCLUSAL</b>	SOFTWARE VERSION	1.5	

### 13. FAULT MESSAGES

The following chart gives a list of fault messages that may appear while the "X-Mind® AC" radiographic system is working.

The chart also includes the causes of the fault messages and how to solve them.

Fault Message	Cause	Solution
E00	RX1 tube head is NOT connected or is out of order	Call the Technical Support
E01	RX2 tube head is NOT connected or is out of order	Call the Technical Support
E02	Corrupted EEPROM data	Call the Technical Support
E03	EEPROM data not saved properly	Call the Technical Support
E05	Line voltage value not included within the set limits	Call the Technical Support
E07	Line voltage value not included within the -5% +10% nominal value.	Call the Technical Support
E08	The "X-RAY" button always seems to be pressed	Make sure it is not jammed
E09	Anomaly in the CONTROL PANEL	Call the Technical Support
E12	The exposure has been prematurely interrupted	Hold the "X-RAY" button down till the end of the exposure
E20	Anomaly in the triac/relay	Call the Technical Support
E21	Anomaly in the electronic circuit	Call the Technical Support
E22	Anomaly in the control circuit	Call the Technical Support
E23	Incorrect dip-switch configuration setting	Call the Technical Support
E24	The "X RAY" button does NOT correspond to the selected tube head	Select the "X-RAY" button that corresponds to the selected tube head or ask Technical Support to verify the configuration.
Err	Indicates a "MAJOR FAULT"	All radiographic system functions are disabled. Call the Technical Support.

## 14. REPLACING FUSES

4 fuses located on the electronic card protect the switchboard electronic equipment.

To replace them proceed as follows (Fig. 25):

1. Isolate the power supply
2. Temporarily remove the guard of the switchboard ① by unscrewing the fixing screws ②


Fig. 25

3. Locate the fuse to be replaced (Fig. 26)
4. Remove the plastic protection (Fig. 27)
5. Remove the fuse
6. Replace it with one of the same type

Power voltage of the radiographic system	220 V -5% +10%	230 V -5% +10%	240 V -5% +10%	115 V -5% +10%
Internal protection fuses (5x20) quick	6.3AF – 250V	6.3AF – 250V	6.3AF – 250V	6.3AF – 250V

7. Refit the cover

Fig. 26


Fig. 27

8. Close the switchboard guard
9. Turn on the power

## **15. CLEANING THE OUTER SURFACES**

Use a soft damp cloth and soap and water to clean the outer surfaces.

The spacer cone may be cleaned with cotton wool soaked in surgical alcohol.

## 16. MAINTENANCE

To ensure the radiographic system works safely, the following checks must be carried out:

### CAUTION

***The owner is responsible for organizing and respecting the maintenance schedule. The radiographic system should be serviced every 12 months to ensure it is working correctly.***

1. Isolate the power
2. Release the spring of the arm "B" of the pantograph using the key provided  
(see §5.7 "BALANCING THE PANTOGRAPH TYPE ARM")
3. Remove the radiographic unit  
(see §5.6 "ASSEMBLING THE RADIOGRAPHIC UNIT")
4. Withdraw the wall plate guard  
(see §5.2 "ASSEMBLING THE WALL PLATE")
5. Remove the terminal board cover and disconnect the pantograph cable
6. Remove the bracket plug and withdraw the guard slab  
(see §5.4 "ASSEMBLING THE PANTOGRAPH TYPE ARM")
7. Remove the pantograph and the relevant cable from the bracket
8. Remove the bracket from the wall plate  
(see §5.3 "ASSEMBLING THE BRACKET")
9. Check the vertical alignment of the wall plate – adjust if required
10. Check the six fixing screws of the wall plate – tighten if required
11. Clean the old lubricating grease from the bracket shaft
12. If the bracket shaft is damaged, install a new bracket
13. Clean the old lubricating grease off the bracket bush
14. If the bracket bush is damaged, install a new bracket
15. Lubricate the bracket shaft (use MOLIKOTE D grease)
16. Lubricate the bracket bush with lubrication grease (use MOLIKOTE D grease)
17. Install the bracket in the wall plate  
(see §5.3 "ASSEMBLING THE BRACKET")
18. Check the pantograph cable – if it is damaged, send the pantograph to the manufacturer for repair
19. Check the pantograph guards
20. Replace the damaged guards
21. Clean the old grease off the shaft  
(see §5.4 "ASSEMBLING THE PANTOGRAPH TYPE ARM")
22. If the shaft is damaged, send the pantograph to the manufacturer for repairs
23. Lubricate the pantograph shaft with lubricating grease (use grease MOLIKOTE D) and reposition it in the bracket
24. Put the pantograph cable back into the bracket and the wall plate, connect it to the terminal board and refit the terminal board cover  
(see §5.5 "CONNECTION TO THE TERMINAL BOARD")
25. Position the guard slab in the bracket
26. Position the bracket plug
27. Position the plate guard
28. Check the electrical contact of the radiographic unit
29. If damaged, send the radiographic unit to the manufacturer for repair
30. Clean the old grease from the radiographic unit assembly shaft
31. Lubricate the assembly shaft of the radiographic unit with a thin layer of lubricating grease (use MOLIKOTE D grease)
32. Reposition the radiographic unit
33. Load the spring of the pantograph arm "B" using the key provided  
(see §5.7 "BALANCING THE PANTOGRAPH TYPE ARM")
34. Power up and check that the radiographic system is working correctly

## 17. REPAIRS

In the event of a breakdown, send the faulty part, (**USING THE ORIGINAL PACKAGING**) to:


SATELEC® S.A.S.  
Z.I. du Phare  
B.P. 216  
33708 MERIGNAC CEDEX  
FRANCE

Tél. + 33 (0) 556 34 06 07  
Fax + 33 (0) 556 34 92 92

E-mail : [info@satelec.com](mailto:info@satelec.com)

## 18. ATTACHMENTS

The manufacturer shall, upon request, supply any drawings, circuit diagrams, component parts lists, instructions or other information needed by qualified technical personnel to carry out repairs on those parts of the "X-Mind® AC" radiographic system which can be repaired.

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